HEIMELIG  
BETTEN

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# Patientenlift Eagle 625

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# VERMEIREN

## 620 625 - Patienthoist

INSTRUCTION MANUAL  
MODE D'EMPLOI  
GEBRUIKSAANWIJZING  
GEBRAUCHSANWEISUNG  
ISTRUZIONI PER L'USO  
MANUAL DE INSTRUCCIONES  
INSTRUKCJA OBSŁUGI

625

620

## **Instructions to the specialist dealer**

This instruction manual is part and parcel of the product and must accompany every product sold.

Version: B, September 2011

All rights reserved, including translation.

No part of this manual may be reproduced in any form what so ever (print, photocopy, microfilm or any other process) without written permission of the publisher, or processed, duplicated or distributed by using electronic systems.

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---

**Contents**

|                                                                                                   |           |
|---------------------------------------------------------------------------------------------------|-----------|
| <b>Preface.....</b>                                                                               | <b>2</b>  |
| <b>1 Product description .....</b>                                                                | <b>3</b>  |
| 1.1 Intended Use.....                                                                             | 3         |
| 1.2 Technical specifications .....                                                                | 4         |
| 1.3 Drawing.....                                                                                  | 6         |
| 1.4 Explanation of the symbols .....                                                              | 6         |
| 1.5 Delivery.....                                                                                 | 7         |
| <b>2 Use.....</b>                                                                                 | <b>7</b>  |
| 2.1 General notes.....                                                                            | 7         |
| 2.2 Use of the battery charger → Eagle 620 .....                                                  | 8         |
| 2.3 Use of the JUMBO system (battery, control box + included battery charger)→<br>Eagle 625 ..... | 9         |
| 2.4 Operating the lift .....                                                                      | 10        |
| 2.5 Slings.....                                                                                   | 12        |
| 2.6 Safety rules.....                                                                             | 12        |
| 2.7 Emergency.....                                                                                | 12        |
| <b>3 Installation and adjustment.....</b>                                                         | <b>14</b> |
| 3.1 Swapping the battery .....                                                                    | 14        |
| 3.2 Assembly or dismantling.....                                                                  | 14        |
| <b>4 Maintenance .....</b>                                                                        | <b>17</b> |
| 4.1 Regular Maintenance.....                                                                      | 17        |
| 4.2 Tools .....                                                                                   | 17        |
| 4.3 Shipping and storage .....                                                                    | 18        |
| 4.4 Care .....                                                                                    | 18        |
| 4.5 Inspection .....                                                                              | 19        |
| 4.6 Troubleshooting.....                                                                          | 20        |
| 4.7 Disinfection .....                                                                            | 21        |
| <b>5 Guarantee .....</b>                                                                          | <b>23</b> |
| <b>6 Disposal.....</b>                                                                            | <b>23</b> |
| <b>7 Declaration of conformity .....</b>                                                          | <b>23</b> |
| <b>8 Maintenance plan .....</b>                                                                   | <b>24</b> |
| <b>9 Disinfection book .....</b>                                                                  | <b>24</b> |
| <b>10 Ordernumbers .....</b>                                                                      | <b>25</b> |

## Preface

First of all we want to thank you for putting your trust in us by selecting one of our products.

The Vermeiren patient hoists are the result of many years of research and experience. During the development, special attention was given to the ease of use and the serviceability of the patient hoist.

The expected lifetime of your patient hoist is strongly influenced by the care and maintenance of the patient hoist.

This manual will help you get acquainted with the operation of your product.

Following the user instructions and the maintenance instructions are an essential part of the guarantee.

This manual reflects the latest product developments. Vermeiren has the right to introduce changes without the obligation to adapt or replace previously delivered models.

For any further questions, please consult your specialist dealer.

## 1 Product description

### 1.1 *Intended Use*

The patient hoist is intended for transferring sitting patients between e.g. wheelchairs, nursing beds and bathrooms.

The patient hoist is designed to transport 1 person.

The patient hoist is suited for indoor use.

The user in the patient hoist must be lifted by an attendant.

The hoists features allow it to be used indoors e.g. for:

- paralysis
- loss of limbs (leg amputation)
- limb defects / malformation
- stiff or damaged joints
- cardiovascular problems
- cachexia
- and also for aged persons.

When providing for individual requirements:

- body size and weight (max. 150 kg for Eagle 620 and 175 kg for Eagle 625)
- physical and psychological condition
- residential circumstances
- environment

should be taken into consideration.

Your patient hoist should only be used on flat surfaces and with all four casters in contact with the ground.

It must not be used in wet rooms/bathrooms or humid environments.

Riding over obstacles while transferring the patient is completely forbidden.

The patient hoist should not be used as a ladder, nor is it a transport for heavy or hot objects.

When used on mats, carpeted floors or loose floor coverings, the floor covering can get damaged.


Use only Vermeiren approved accessories.

Make sure that the end of the adjustment motor has a safety measurement. Therefore, use the patient hoist only in the middle of the adjustment range and not at the end of the hoist motor.

Compliance with the user and maintenance instructions are an essential part of the guarantee conditions.

Visually impaired people can contact the dealer for the instructions for use.

## 1.2 Technical specifications

| Brand                                              | Vermeiren                             |                                                         |                          |
|----------------------------------------------------|---------------------------------------|---------------------------------------------------------|--------------------------|
| Address                                            | Vermeirenplein 1/15, B-2920 Kalmthout |                                                         |                          |
| Type                                               | Patient hoist                         |                                                         |                          |
| Model                                              | Eagle                                 |                                                         |                          |
|                                                    | Measurement in drawing                | Dimensions 620                                          | Dimensions 625           |
| Lowest CAP position*                               | l                                     | 728 mm                                                  | 766 mm                   |
| Highest CAP position*                              | k                                     | 1740 mm                                                 | 1840 mm                  |
| Lifting range (range of heights)                   | m                                     | 1012 mm                                                 | 1074 mm                  |
| Leg length                                         |                                       | 1000 mm                                                 | 1150 mm                  |
| Total length                                       |                                       | 1115 mm                                                 | 1290 mm                  |
| Minimum leg separation                             | r                                     | 463 mm                                                  | 554 mm                   |
| Maximum leg separation                             | q                                     | 833 mm                                                  | 945 mm                   |
| Leg height / Frame height                          | n                                     | 111.5 mm                                                | 111.5 mm                 |
| Overall width (closed), external dimension         |                                       | 550 mm                                                  | 680 mm                   |
| Overall width (opened), castors forwards           | p                                     | 1035 mm                                                 | 1225 mm                  |
| Min. distance between wall / CSP* (minimum height) | w3                                    | 680 mm                                                  | 750 mm                   |
| Min. distance between wall / CSP* (maximum reach)  | w2                                    | 240 mm                                                  | 310 mm                   |
| Min. distance between wall / CSP* (maximum height) | w1                                    | 320 mm                                                  | 390 mm                   |
| Turning circle                                     |                                       | 1160 mm                                                 | 1290 mm                  |
| Total weight                                       |                                       | 37 kg                                                   | 42 kg                    |
| Weight of chassis + motor                          |                                       | 17 kg                                                   | 20 kg                    |
| Weight of battery                                  |                                       | Not applicable                                          | 2.90 kg                  |
| Weight of cantilever arm / expanding bracket       |                                       | 20 kg                                                   | 22 kg                    |
| <b>Maximum load</b>                                |                                       | <b>150 kg</b>                                           | <b>175 kg</b>            |
| Free height, minimum                               |                                       | 51.5 mm                                                 | 51.5 mm                  |
| Maximum reach at 600 mm                            | a                                     | 265 mm                                                  | 434 mm                   |
| Maximum reach from support                         | b                                     | 476 mm                                                  | 508 mm                   |
| Reach from support at separation 700 mm            | c                                     | 295 mm                                                  | 172 mm                   |
| Power outlet                                       |                                       | 24V --- max. 250 VA                                     |                          |
| Supply voltage                                     |                                       | 24V --- max. 5.3 VA                                     | 100-240V ~ max. 37-53 VA |
| Maximum current drawn                              |                                       | max. 300 mA                                             | max. 400 mA              |
| Operating temperature                              |                                       | +5 to +40°C                                             |                          |
| Relative humidity                                  |                                       | 20% to 90% at 30°C, non-condensing                      |                          |
| Noise level                                        |                                       | < 50 dB(A)                                              |                          |
| Air pressure                                       |                                       | 700 to 1060 hPa                                         |                          |
| Control box                                        |                                       | Linak CBJ Home                                          | Linak CBJ2               |
| Battery                                            |                                       | Not applicable                                          | Linak BAJ1               |
| Manual control                                     |                                       | Linak HB7X (< 5N force required to operate)             |                          |
| Motor                                              |                                       | Linak LA31C (6000N)                                     | Linak LA34 (10000N)      |
| Control box protection class                       |                                       | IPX4                                                    |                          |
| Battery protection class                           |                                       | Not applicable                                          | IPX5                     |
| Manual control protection class                    |                                       | IP66                                                    |                          |
| Motor protection class                             |                                       | IP54                                                    |                          |
| Insulation class                                   |                                       | II - Type B                                             |                          |
| Operational performance                            |                                       | approx. 40 lifts per charge                             |                          |
| Periodic operation                                 |                                       | max. 10%, or 2 min. continuous operation/18 min. pauses |                          |
| Battery capacity                                   |                                       | 2.9 Ah                                                  |                          |
| Emergency stop                                     |                                       | Yes                                                     |                          |

|                                  |                                       |                |                |
|----------------------------------|---------------------------------------|----------------|----------------|
| Brand                            | Vermeiren                             |                |                |
| Address                          | Vermeirenplein 1/15, B-2920 Kalmthout |                |                |
| Type                             | Patient hoist                         |                |                |
| Model                            | Eagle                                 |                |                |
|                                  | Measurement in drawing                | Dimensions 620 | Dimensions 625 |
| Manual lowering in emergency     |                                       | No             | Yes            |
| Electrical lowering in emergency |                                       | Yes            | No             |

We reserve the right to introduce technical changes. Measurement tolerance  $\pm 15$  mm / 1,5 kg

1 = Highest position, 2 = Maximum reach, 3 = Lowest position


\* CSP = Central suspension point

Table 1: Technical specifications

The patient hoist complies to the requirements set up in:

EN ISO 10535: Hoists for the transfer of disabled persons - Requirements and test methods

### 1.3 Drawing

### 1.4 Explanation of the symbols

Observe the safety instructions!

Read the instruction manual before use!

Separate recovery and recycling of electric and electronic devices.

CE declaration of conformity

## 1.5 Delivery

The Vermeiren Eagle 620 patient hoist shall be delivered with:

- chassis including 4 casters (2 with brakes)
- pole, including handles for pushing
- cantilever arm and expanding bracket
- control box (incl. 2 accumulator batteries and manual operation)
- battery charger
- motor
- instruction manual

The Vermeiren Eagle 625 patient hoist shall be delivered with:

- chassis including 4 casters (2 with brakes)
- pole, including handles for pushing
- cantilever arm and expanding bracket
- control box (+ included battery charger)
- battery
- manual control
- motor
- instruction manual

Before use check if everything is included and that no products are damaged (example by transport, ...).

Please note that the basic configuration may vary from one European country to another. Please contact the specialist dealer in your country for more information.

## 2 Use

This chapter describes the everyday use. **These instructions are for the user and the specialist dealer.**

The patient hoist is delivered fully assembled by your specialist dealer. The instructions intended for the specialist dealer on how to set up the patient hoist are given in § 3.

### 2.1 General notes

The patient hoist must only be operated by qualified staff who have been instructed or trained in its specific application.

The use of the patient hoist is restricted to indoors only. Please ensure that it is only used on flat surfaces. It must be used strictly indoors only. Use only the battery charger (Eagle 620) supplied and no other charging equipment.

We would like to point out that sources of electromagnetic radiation (e.g. mobile phones, etc.) can cause interference and that the lift's electronics can also affect other electrical appliances.

Even if you have been instructed by your dealer about the operational elements of your lift and their use, we recommend that you read the following pages carefully.

## 2.2 Use of the battery charger → Eagle 620

### 2.2.1 Battery charger

|                                      |                                                                                   |
|--------------------------------------|-----------------------------------------------------------------------------------|
| Primary voltage                      | 100 - 240VAC / 50/60 Hz                                                           |
| Secondary voltage (charging voltage) | 27.6VDC (+/- 2%)                                                                  |
| Secondary current (charging current) | max. 500 mA                                                                       |
| Safeguards                           | Protected against reverse polarisation, electrical surges and extreme temperature |
| Efficiency                           | min. 80% (when fully loaded)                                                      |
| Ambient temperature                  | 0 to +40°C                                                                        |
| Length of mains lead                 | 2.0 m                                                                             |
| Ambient storage temperature          | -15° to +50°C                                                                     |
| Relative storage humidity            | max. 95% (non-condensing)                                                         |
| Conformity                           | tested in accordance with IEC 60601-1                                             |

### 2.2.2 Charging the batteries

Use only the battery charger supplied and no other charging equipment.

We recommend that you should regularly recharge the built-in control box batteries in order to ensure that the hoist remains usable and to extend the lifespan of the batteries. The control box will give the warning sound to tell you when the battery charge is too low.

- **SETTING UP THE CHARGER**

When setting up the charger, ensure that it is well-ventilated on all sides. A minimum of 100 mm space should be left free around the unit for this purpose. If the charger is insufficiently ventilated so that the unit heats up, the charging rate will fall which will extend the charging time. If the charger overheats ( $> +50^{\circ}\text{C}$ ), it will stop charging.

The charger should only be used from a wall socket with a voltage of 100 to 240V at 50/60Hz and in well-ventilated, dry areas.

- **FIRST USE**

First, insert the recharge cable into the appropriate connector on the control box. Then insert the mains adapter into the socket. The loading time is about 4 hours. The charger will stop the recharging process automatically when the device is fully loaded.

- **RECHARGING**

 **WARNING: Risk of injury – Do not use the lift when the charger is connected in the control box.**

When the batteries are recharging there will burn a yellow light on the upper LED.

After the recharging process is completed, you should always remove the mains plug of the charger from the socket first and then remove the charger from the control box.

Do not use the lift while it is being recharged!

If the batteries are not used for a significant period, they discharge slowly by themselves (deep discharge). It then becomes impossible to recharge them with the battery charger supplied. You should therefore recharge the batteries at least once a month even if they are not being used.

## 2.3 Use of the JUMBO system (battery, control box + included battery charger)→ Eagle 625

JUMBO is a modular system combining an actuator, a control box and a battery in a flexible solution, specially developed for patient hoists.

### 2.3.1 Battery

The battery is above the control box and is replaceable with an integrated clip system.

|                             |                                       |
|-----------------------------|---------------------------------------|
| Primary voltage             | 24 VDC                                |
| Ambient temperature         | +5°C to +40°C                         |
| Ambient storage temperature | -10° to +50°C                         |
| Relative storage humidity   | max. 90% (non-condensing)             |
| Conformity                  | tested in accordance with IEC 60601-1 |

### 2.3.2 Control box + included battery charger

The control box is designed with a red emergency stop and an internal battery charger.

|                                      |                                                                                   |
|--------------------------------------|-----------------------------------------------------------------------------------|
| Primary voltage                      | 100 - 240 VAC / 50/60 Hz                                                          |
| Secondary voltage (charging voltage) | 24 VDC, max. 250 VA                                                               |
| Secondary current (charging current) | max. 10 A                                                                         |
| Safeguards                           | Protected against reverse polarisation, electrical surges and extreme temperature |
| Ambient temperature                  | +5°C to +40°C                                                                     |
| Ambient storage temperature          | -10° to +50°C                                                                     |
| Relative storage humidity            | max. 90% (non-condensing)                                                         |
| Conformity                           | tested in accordance with IEC 60601-1                                             |

We reserve the right to introduce technical changes.

### 2.3.3 External charger (OPTIONAL)

The external charger shall be mounted on the wall. It is possible to charge the spare battery (optional) with this external charger and you have no waiting time to use the patient hoist when charging the batteries.

|                             |                           |
|-----------------------------|---------------------------|
| Primary voltage             | 100 - 240 VAC / 50/60 Hz  |
| Charging current            | max. 650 mA               |
| Ambient temperature         | +5°C to +40°C             |
| Ambient storage temperature | -10° to +50°C             |
| Relative storage humidity   | max. 90% (non-condensing) |

### 2.3.4 Charging the batteries

Use only the control box with included battery charger on the Eagle 625 or the external charger meant to be mounted on the wall.

We recommend that you should regularly recharge the batteries in order to ensure that the hoist remains usable and to extend the lifespan of the batteries. The control box will give a warning sound to tell you when the battery charge is on 50% capacity.

- **FIRST USE**

First, insert the supply cable with the connector side into the appropriate connector on the control box. The plug side of the supply cable must be insert in the socket outlet. The loading time is about 24 hours.

- **RECHARGING**

**⚠ WARNING: Risk of injury – Do not use the lift when the supply cable is connected in the socket outlet.**

When the batteries are recharging there will burn a yellow light (above on) and an orange light (above charge).

After the recharging process is completed, you should always remove the mains plug from the socket first and then remove the connector plug from the control box.

Do not use the lift while it is being recharged!


If the batteries are not used for a significant period, they discharge slowly by themselves (deep discharge). It then becomes impossible to recharge them. You should therefore recharge the batteries at least once a month even if they are not being used.

## 2.4 Operating the lift

Please observe the technical information according to which the patient hoists (620 or 625) may be operated. The patient hoist must only be operated by authorised staff who have been trained in its use and operation.

**⚠ WARNING: Risk of burns – Be careful when driving in hot or cold environments (sunshine, extreme cold, saunas, etc.) for a sufficient amount of time and when touching - Surfaces can assume environment temperatures.**


### 2.4.1 Altering the chassis width

To get the lift around a wheelchair or other seating furniture, or to increase the stability of the lift in the standing position, the separation of the chassis legs can be increased.


Stand behind the assembled lift and take hold of the handles for pushing the patient hoist (on the left and right, next to the control box). Press the lever at the bottom of the chassis (left or right) down with your foot, and the separation of the legs of the chassis can be increased or decreased.

## 2.4.2 Parking brakes

Secure the two casters at the back of the chassis by pressing the caster brake plates down gently to their end stops with the tip of your foot. To release the brake, push the brake plate gently back up again with the front of your foot until the casters are free.

## 2.4.3 Raising and lowering the cantilever arm

The manual control allow the cantilever arm to be adjusted smoothly to any point.

- ① = Raise the cantilever arm
- ② = Lower the cantilever arm

⚠ When using the lift function until the end of the hoist motor the limit switch is activated by electronics (security function). Use the lifting function only in the middle of the adjustment range.

Force required to activate the pushbutton: < 5N

Position the manual control to the top tube of the handles for moving the patient hoist. This is the most ergonomic position.

## 2.4.4 Moving the patient hoist

⚠ **WARNING: Risk of clamping – Be careful passing through restricted passages (e.g. doors).**

1. Stand behind the patient hoist and grasp the handles (left and right, next to the control box) with both hands.
2. Make sure that both back caster brakes have been released.
3. Push the patient hoist slowly into the desired position.

When it is not loaded, we advise you to pull the patient hoist backwards, which makes it easier to get round any obstacles (door frames, corners of rooms or furniture).

**To move the patient hoist use only the handles for pushing and no other components (actuator, manual control, ...) .**

## 2.5 Slings

- ⚠ **WARNING: Risk of injury – Only use appropriate slings for the patients.**
- ⚠ **WARNING: Risk of injury – Do not use damaged slings.**

The Eagle patient hoist must only be used with slings designed for the Vermeiren Eagle patient hoists (for 4-point expanding brackets).

The applicability and instructions for use for the specific slings should be followed.

## 2.6 Safety rules


Some safety tips are given below for your security:

- ⚠ Only authorised staff who have been trained in using the lift must use and operate the patient hoist.
- ⚠ Make sure when using the lift that there is sufficient room around and above it, as the adjustment movements could otherwise lead to damage or injury.
- ⚠ Always be aware of the lifting arm to avoid injury.
- ⚠ The patient hoist must only be used on level surfaces, with all 4 casters contacting the ground evenly.
- ⚠ It must not be used in humid environments.
- ⚠ Only slings that have been designed and approved for the patients must be used (see usage instructions for the various slings). Use of any other slings is at your own risk.
- ⚠ Please follow the instructions of the nursing staff or suitably trained persons, so that no injuries are caused during use of the lift.
- ⚠ Check the medical condition of the patient and be sure that it is possible to lift the patient with this patient hoist. (Refer to intended use)
- ⚠ Be careful with sources of ignition such as lit cigarettes as they may set the sling alight.
- ⚠ Never exceed the maximum load (**150 kg for the 620 and 175 kg for the 625**). The control box is disabled if the lift is overloaded.
- ⚠ When using slings with a maximum load, the lowest load is used between the patients lift and sling.
- ⚠ The manufacturer accepts no liability for damages and injury caused by inappropriate actions or by not observing the instructions given in this manual.
- ⚠ Use only the handles for pushing to move the patient hoist and no other components (actuator, manual control, ...).

## 2.7 Emergency

In cases for emergency the emergency stop and the function for lowering in an emergency must be used.

### 2.7.1 Emergency stop (Eagle 620)

In emergencies, and during assembly and disassembly, the emergency stop is activated by pressing the red button on the right of the control box. Push on the button and turn in the direction of the arrows. If the button is fully pressed to the back turn back in the opposite direction of the arrows.


The emergency stop is deactivated by turning the red button in the direction of the arrows.

## 2.7.2 Emergency stop (Eagle 625)

In emergencies, and during assembly and disassembly, the emergency stop is activated by pressing the red button on the control box.

Push on the button and turn in the direction of the arrows. If the button is fully pressed to the back turn back in the opposite direction of the arrows.

The emergency stop is deactivated by turning the red button in the direction of the arrows.

## 2.7.3 Lowering in an emergency (Eagle 620)


If the manual operation mechanism should be damaged, it is possible to lower the cantilever arm from the control box.

The arm moves down for as long as the bottom button on the control box panel is pressed.

## 2.7.4 Manual emergency lowering (Eagle 625)

**⚠ WARNING: Risk of injury – Adjust the manual emergency lowering for the weight of the patient.**

Emergency lowering if the power fails or the batteries are empty can be done using the red pull button ① at the lower end of the motor. This emergency lowering is factory-set for a patient weight of 75 kg.

Standard calibration: 3000N, 16 mm/s to lower the patient.

Please note that the manual emergency lowering is only possible when the patient is sitting in the patient hoist.

Handle pull force: 10 mm for force from 80 à 100N.

### 3 Installation and adjustment

The instructions in this chapter are for the specialist dealer.


To find a service facility or specialist dealer near you, contact the nearest Vermeiren facility. A list of Vermeiren facilities can be found on the last page.

**⚠ WARNING: Risk of unsafe limitations - Use only the limitations described in this manual.**

#### 3.1 Swapping the battery

- We shall decline all liability for damage caused by the use of improper battery.
- Do not use the battery at temperatures below +5°C or above +50°C (the ideal is +20°C).
- If the battery is opened, all liability of the manufacturer shall be voided as well as any claim.

The battery of the patient hoist Eagle 625 can be swapped by an integrated clip system.


To swap the battery pull on the integrated clip system. The plate will be lifted and the battery can take off the holder (control box).

#### 3.2 Assembly or dismantling


##### 3.2.1 Unpacking

**⚠ WARNING: Risk of injury - After unpacking and before further assembly, always first confirm that the emergency stop (the red pushbutton on the control box) is pressed in.**

The patient hoist packaging has been selected to provide the best possible protection during transport.


1. Take the patient hoist out of the box and check that all the items have been delivered and that the individual parts do not have any visible defects. If there is any damage, please contact your sales office.


2. Before assembly, put on the caster brakes to make sure that the lift cannot be moved unintentionally. Push the caster brake plate down gently to its end stop with the front of your foot. To release the brake, push the brake plate gently back up again with the front of your foot until the caster is free.

### 3.2.2 Assembly


**⚠ WARNING: Risk of clamping - Take care that no body parts and wires get trapped, crushed or cut during assembly.**


1. Remove the safety spindle (locking pin) from the bottom end of the pole by gently pressing the button at the head of the spindle. The safety spindle can now be removed easily.

2. Undo the transport lock (wing screw) gently and the pole can now be raised.


- 2.a The standing position safety (traction screw) clicks audibly into place once the pole is put in its final position.

3. **⚠ WARNING: Risk of injury - Make sure that the safety spindle is correctly inserted.**


Fix the pole in place by pushing the safety spindle at the end of the pole in with the pushbutton depressed until its other end clearly protrudes. Release the pushbutton on the safety spindle and check that it cannot be pulled out.


4. Undo the transport strap (Velcro) of the expanding bracket. This can now be moved freely to the end of the pole.

### 3.2.3 Dismantling


- ⚠ **WARNING: Risk of clamping - Take care that no body parts and wires get trapped, crushed or cut during dismantling.**
- ⚠ **WARNING: Risk of injury - Check the emergency stop (red pushbutton on the control box) is activated, to avoid unintentional adjustment movements.**


1. Secure the two casters at the back of the chassis by pressing the caster brake plates down gently to their end stops with the tip of your foot. To release the brake, push the brake plate gently back up again with the front of your foot until the casters are free.

2. Move the cantilever arm to the bottom position (end stop). Place the expanding bracket vertically and secure it with the transport strap Velcro, as shown in the picture.


3. Remove the safety spindle (locking pin) from the bottom end of the pole by gently pressing the button at the head of the spindle. The safety spindle can now be removed easily.

4. ⚠ **WARNING: Risk of injury - Once the pole is unlocked, its own weight can cause considerable downward forces to be exerted.**

Pull out the standing position safety (traction screw) until the pole is released and can be tilted. The pole can now be tilted forwards.

Tip the pole as far as possible forward that its head end comes to rest on the expanding bracket.

5. ⚠ **WARNING: Risk of injury - Make sure that the safety spindle is correctly inserted.**

Fix the pole in place by pushing the safety spindle at the end of the pole in with the pushbutton depressed until its other end clearly protrudes. Release the pushbutton on the safety spindle and check that it cannot be pulled out again.

## 4 Maintenance

The expected lifetime of the patient hoist is 8 years, depending on type of use, storage, regular maintenance, servicing and cleaning.

### 4.1 Regular Maintenance

The maintenance activities to ensure a good condition for your Vermeiren patient hoist are described below:

#### ● BEFORE EACH USE

1. Check the brakes for visible damage and/or dirt. Remove any dirt as it could impair the motion and grip of the wheels. If a caster is damaged, please go to an authorised workshop to get it repaired.
2. Check that the electrical adjustment of the cantilever arm is working properly, using the manual control. If this functionality is impaired, contact your specialist dealer.
3. Ensure that the batteries are sufficiently loaded for the intended use and recharge them as required.
4. Check the patient hoist (cleanliness, cracks, damage on structural parts, ...) and make it clean. If necessary change the coating.

#### ● APPROX. EVERY 8 WEEKS

Depending on the frequency of use, check the following:

1. Lubrication of the joints of the cantilever arm
2. Condition of the casters
3. Whether there is any visible damage to the casing of the control box, battery, battery charger and the manual controls
4. Whether the cables are fit for use
5. Recharge the battery

#### ● APPROX. EVERY 6 MONTHS OR FOR EVERY NEW USER

Depending on the frequency of use check the following:

1. Cleanliness
2. General condition
3. That the battery charger is working (Eagle 620)
4. Casters are working

If the resistance to rotation is too great, clean the casters. If this is insufficient, please consult your dealer.

For your convenience, on the back of this manual is a maintenance plan.

Repairing and assembling of spare parts for your patient hoist may only be performed by the specialist dealer.

Only authorized Vermeiren spare parts may be installed.

### 4.2 Tools

**⚠ CAUTION: Risk of injury - Any activities requiring tools should only be carried out by authorized persons.**

The following tools are the minimum required for repairing individual components and accessory attachments, and for carrying out inspections:

- Allen keys (sizes 3, 4, 5 and 6 mm)
- Fixed spanners/ring spanners (sizes 8, 10, 12, 13, 17 and 19 mm)
- Phillips screwdrivers (PH1, PH2)

### 4.3 Shipping and storage

The shipping and storage of the patient hoist shall be according following instructions:

- Store in a dry place (between + 5°C and + 41°C).
- The relative humidity should be between 30% and 70%.
- Check cables are not being crushed and to prevent kinks.
- Disconnect the power plug of the battery charger (Eagle 620).
- Disconnect the patient hoist out the socket outlet (Eagle 625).
- Provide sufficient covering or packaging to protect the patient hoist from rust and foreign bodies. (e.g. salt water, sea air, sand, dust).
- Store all removed parts together in one place (or mark them if necessary) to avoid mixing up with other products when re-assembling (e.g. the charger).
- Components must be stored without being subjected to strains (Do not put too heavy parts on the patient hoist, not clamping between something, ...).

### 4.4 Care

Your lift will also require regular care if you want it to remain visually attractive. Please take note of the following tips:

- Use a mild commercial detergent for removing stubborn dirt.
- Stains can be removed by using a sponge or a soft brush.
- Do not use strong cleaning liquids like solvents, nor use hard brushes.
- Never clean with steam and/or pressure washers.

#### 4.4.1 Plastic parts

Clean plastic parts of your patient hoist with commercial plastic cleaners. Read the specific product information and only use a soft brush or soft sponge.

#### 4.4.2 Coating

The high quality of the surface layer guarantees optimal protection against corrosion. If the outer coat is damaged by scratching or in some other way, get your specialist dealer to repair the affected surface.

When cleaning, only use warm water and normal household detergents and soft brushes and cloths. Ensure that no wetness gets into the tubes.

#### 4.4.3 Electronics casing

**⚠ WARNING: Risk of injury - The emergency stop must be activated before maintenance work, so that no unintentional adjustment movements can be made.**

Wipe the control box, motor casing and manual controls down with a cloth moistened by a few drops of a commercial domestic cleaner. Do not use abrasives or sharp-edged polishing equipment like a metal scrubber or brush, as these can scratch the surface of the manual control and destroy its water repelling property.

Regularly check the plug connectors for corrosion or damage, as these could affect the functional integrity of the electronics.

The manufacturer will not accept liability for damages caused by insufficient maintenance.

## 4.5 Inspection

In principle we recommend one inspection every year, and a minimum of one before usage is resumed. All of the following checks must be performed and documented by authorized persons:

- Check the cabling (in particular for crushing, wear, cuts, visible parts of the insulation of the inner wires, visible metal wire cores, kinks, bulges, colour changes of the outer sleeve, brittle spots, and safe positioning so that mechanical effects such as being cut or crushed are not probable).
- Visual inspection of the frame parts to check for plastic deformation and/or wear and tear (frame, motor suspension, cantilever arm, expanding bracket).
- Visually check for damage to the paintwork (danger of corrosion).
- Visual inspection of all casings for damage; screws must be tight.
- Verify the amount of grease on the metal joints of movable parts
- Visually check all plastic parts for cracks and brittle spots
- Review the charger (Eagle 620), control box + included charger (Eagle 625) residual discharge current (A) based on VDE 0702
- Review the charger (Eagle 620), control box + included charger (Eagle 625) insulation resistance (MO) based on VDE 0702
- Functional test of the cantilever arm (lubrication of the joint parts, range of adjustability, deformation, wear and tear)
- Check the functioning of the hoist drives (during a test drive → noises, speed, free running, etc.), if necessary: Measure the performance, first with no load and then with the nominal load ("SWL"), to investigate the wear and tear of the motors by comparing the values for the electric current against the values when it was delivered.
- Functional test of the emergency stop
- Functional test of the emergency lowering
- Completeness of the delivery condition, instruction manual available?

Checking measurements may only be carried out by skilled persons trained on the patient hoist at least and at least under the supervision of an electrician who knows the checking instruments and processes. Only an electrician can release the patient hoist for use after making the measurements and carrying out the servicing.

The service must only be signed off in the maintenance plan if a minimum of all the above-mentioned aspects have been checked.

## 4.6 Troubleshooting

| Problem                                                              | Cause of problem                                                                                                                                                                                                                                                                                                                                    | Problem solving                                                                                                                                                                                                                                                                                                                                                                                          |
|----------------------------------------------------------------------|-----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Noisy sound of the moving parts (example: cantilever arm).           | Lubrication is necessary.                                                                                                                                                                                                                                                                                                                           | Lubricate the moving parts.<br><b>(Do not lubricate the actuator!)</b>                                                                                                                                                                                                                                                                                                                                   |
| Actuator makes an unusual sound.                                     | The actuator is damaged.                                                                                                                                                                                                                                                                                                                            | Replace the actuator or contact your specialist dealer.                                                                                                                                                                                                                                                                                                                                                  |
| Patient hoist does not lift anymore.                                 | <ol style="list-style-type: none"> <li>1. Electric actuator is damaged.</li> <li>2. The user weight is too high.</li> <li>3. Manual control or actuator are not connected.</li> <li>4. Battery not connected or not connected properly to control box. (Eagle 625)</li> <li>5. Battery too low.</li> <li>6. Emergency stop is activated.</li> </ol> | <ol style="list-style-type: none"> <li>1. Replace the actuator.</li> <li>2. Use a patient hoist with higher maximum user weight.</li> <li>3. Check connections. If necessary connect manual control or actuator.</li> <li>4. Check battery is installed or the connection of the battery. (Eagle 625)</li> <li>5. Charge the batteries or replace.</li> <li>6. Deactivate the emergency stop.</li> </ol> |
| Lowering the patient in case of emergency does not work. (Eagle 625) | There is required a minimum weight load to lower the patient.                                                                                                                                                                                                                                                                                       | If the patient weight is too low pull down slightly on the cantilever arms or adjust the lowering speed.                                                                                                                                                                                                                                                                                                 |
| Battery cannot be charged                                            | <ol style="list-style-type: none"> <li>1. Batteries not fitted correctly to control box.</li> <li>2. Mains cable defect.</li> <li>3. Batteries are defected.</li> <li>4. Control box with included charger is defected (Eagle 625)</li> <li>5. Battery charger is defected (Eagle 620)</li> </ol>                                                   | <ol style="list-style-type: none"> <li>1. Check batteries are connected to the control box.</li> <li>2. Change mains cable.</li> <li>3. Replace battery.</li> <li>4. Replace control box with included charger.</li> <li>5. Replace new battery charger. (Eagle 620)</li> </ol>                                                                                                                          |

Table 2: Troubleshooting

## 4.7 Disinfection

- ⚠ WARNING: Dangerous Products - The use of disinfectants is restricted to authorized personnel.**
- ⚠ WARNING: Dangerous Products, change on irritate your skin - You should wear suitable protective clothing because the disinfectants could irritate your skin. For this purpose you should also take note of the product information of the solutions concerned.**

All parts of the patient hoist can be treated by scrubbing with a disinfectant.

All steps taken to disinfect rehabilitation equipment, their components or other accessory parts are to be recorded in a disinfection report containing a minimum of the following information (with product documentation appended):

| Date of the disinfection | Reason | Specification | Substance and concentration | Signature |
|--------------------------|--------|---------------|-----------------------------|-----------|
|--------------------------|--------|---------------|-----------------------------|-----------|

Table 3: Example of a disinfection book

Abbreviations used in column 2 (reason):

V = Suspected infection IF = Infection case W = Repetition I = Inspection

For a blank sheet of a disinfection book see § 9.

The recommended disinfectants for scrubbing (based on the list provided by the Robert Koch Institute, RKI) are standing in the table below. The current state of the disinfectants included in the RKI list can be obtained from the Robert Koch Institute (RKI) (homepage: [www.rki.de](http://www.rki.de)).

| Active substance                                               | Product name                     | Laundry disinfection |                     | Surface disinfection (scrubbing-/wiping disinfection) |                     | Disinfection of excretions<br>1 part sputum or stools + 2 parts diluted solution or 1 part urine + 1 part diluted solution |                     |                  |                     |                  |                     | Area of effectiveness | Manufacturer or Supplier |
|----------------------------------------------------------------|----------------------------------|----------------------|---------------------|-------------------------------------------------------|---------------------|----------------------------------------------------------------------------------------------------------------------------|---------------------|------------------|---------------------|------------------|---------------------|-----------------------|--------------------------|
|                                                                |                                  |                      |                     |                                                       |                     | Sputum                                                                                                                     |                     | Stools           |                     | Urine            |                     |                       |                          |
|                                                                |                                  | Diluted solution     | Time to take effect | Diluted solution                                      | Time to take effect | Diluted solution                                                                                                           | Time to take effect | Diluted solution | Time to take effect | Diluted solution | Time to take effect |                       |                          |
|                                                                |                                  | %                    | Hr.                 | %                                                     | Hr.                 | %                                                                                                                          | Hr.                 | %                | Hr.                 | %                | Hr.                 |                       |                          |
| Phenol or phenol derivative                                    | Amocid                           | 1                    | 12                  | 5                                                     | 6                   | 5                                                                                                                          | 4                   | 5                | 6                   | 5                | 2                   | A                     | Lysoform                 |
|                                                                | Gevisol                          | 0,5                  | 12                  | 5                                                     | 4                   | 5                                                                                                                          | 4                   | 5                | 6                   | 5                | 2                   | A                     | Schülke & Mayr           |
|                                                                | Helipur                          |                      |                     | 6                                                     | 4                   | 6                                                                                                                          | 4                   | 6                | 6                   | 6                | 2                   | A                     | B. Braun                 |
|                                                                | m-cresylic soap solution (DAB 6) | 1                    | 12                  | 5                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | A                     |                          |
|                                                                | Phenol                           | 1                    | 12                  | 3                                                     | 2                   |                                                                                                                            |                     |                  |                     |                  |                     | A                     |                          |
| Chlorine, organic or inorganic substances with active chlorine | Chloramin-T DAB 9                | 1,5                  | 12                  | 2,5                                                   | 2                   | 5                                                                                                                          | 4                   |                  |                     |                  |                     | A <sup>1</sup> B      |                          |
|                                                                | Clorina                          | 1, 5                 | 12                  | 2,5                                                   | 2                   | 5                                                                                                                          | 4                   |                  |                     |                  |                     | A <sup>1</sup> B      | Lysoform                 |
|                                                                | Trichlorol                       | 2                    | 12                  | 3                                                     | 2                   | 6                                                                                                                          | 4                   |                  |                     |                  |                     | A <sup>1</sup> B      | Lysoform                 |
| Per combinations                                               | Apesin AP100 <sup>2</sup>        |                      |                     | 4                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Tana PROFESSIONAL        |
|                                                                | Dismozon pur <sup>2</sup>        |                      |                     | 4                                                     | 1                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Bode Chemie              |
|                                                                | Perform <sup>2</sup>             |                      |                     | 3                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Schülke & Mayr           |
|                                                                | Wofesteril <sup>2</sup>          |                      |                     | 2                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Kesla Pharma             |

| Active substance                                   | Product name                                | Laundry disinfection |                     | Surface disinfection (scrubbing-/wiping disinfection) |                     | Disinfection of excretions<br>1 part sputum or stools + 2 parts diluted solution or 1 part urine + 1 part diluted solution |                     |                  |                     |                  |                     | Area of effectiveness | Manufacturer or Supplier |
|----------------------------------------------------|---------------------------------------------|----------------------|---------------------|-------------------------------------------------------|---------------------|----------------------------------------------------------------------------------------------------------------------------|---------------------|------------------|---------------------|------------------|---------------------|-----------------------|--------------------------|
|                                                    |                                             |                      |                     |                                                       |                     | Sputum                                                                                                                     |                     | Stools           |                     | Urine            |                     |                       |                          |
|                                                    |                                             | Diluted solution     | Time to take effect | Diluted solution                                      | Time to take effect | Diluted solution                                                                                                           | Time to take effect | Diluted solution | Time to take effect | Diluted solution | Time to take effect |                       |                          |
|                                                    |                                             |                      |                     |                                                       |                     |                                                                                                                            |                     |                  |                     |                  |                     |                       |                          |
| Formaldehyde and/or other aldehydes or derivatives | Aldasan 2000                                |                      |                     | 4                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Lysoform                 |
|                                                    | Antifect FD 10                              |                      |                     | 3                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Schülke & Mayr           |
|                                                    | Antiseptica surface disinfection 7          |                      |                     | 3                                                     | 6                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Antiseptica              |
|                                                    | Apesin AP30                                 |                      |                     | 5                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | A                     | Tana PROFESSIONAL        |
|                                                    | Bacillocid special                          |                      |                     | 6                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Bode Chemie              |
|                                                    | Buraton 10F                                 |                      |                     | 3                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Schülke & Mayr           |
|                                                    | Desomed A 2000                              |                      |                     | 3                                                     | 6                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Desomed                  |
|                                                    | Hospital disinfectant cleaner               |                      |                     | 8                                                     | 6                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Dreiturm                 |
|                                                    | Desomed Perfekt                             |                      |                     | 7                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB*                   | Desomed                  |
|                                                    | Formaldehyde-solution (DAB 10), (formaline) | 1,5                  | 12                  | 3                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    |                          |
|                                                    | Incidin Perfekt                             | 1                    | 12                  | 3                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Ecolab                   |
|                                                    | Incidin Plus                                |                      |                     | 8                                                     | 6                   |                                                                                                                            |                     |                  |                     |                  |                     | A                     | Ecolab                   |
|                                                    | Kohrsolin                                   | 2                    | 12                  | 3                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Bode Chemie              |
|                                                    | Lysoform                                    | 4                    | 12                  | 5                                                     | 6                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Lysoform                 |
|                                                    | Lysoformin                                  | 3                    | 12                  | 5                                                     | 6                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Lysoform                 |
|                                                    | Lysoformin 2000                             |                      |                     | 4                                                     | 6                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Lysoform                 |
|                                                    | Melsept                                     | 2                    | 12                  | 4                                                     | 6                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | B. Braun                 |
|                                                    | Melsitt                                     | 4                    | 12                  | 10                                                    | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | B. Braun                 |
|                                                    | Minutil                                     | 2                    | 12                  | 6                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Ecolab                   |
|                                                    | Multidor                                    |                      |                     | 3                                                     | 6                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Ecolab                   |
|                                                    | Nüscosept                                   |                      |                     | 5                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB                    | Dr. Nüsken Chemie        |
|                                                    | Optisept                                    |                      |                     | 7                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB*                   | Dr. Schumacher           |
|                                                    | Pursept-FD                                  |                      |                     | 7                                                     | 4                   |                                                                                                                            |                     |                  |                     |                  |                     | AB*                   | Merz                     |
| Ultrasol F                                         | 3                                           | 12                   | 5                   | 4                                                     |                     |                                                                                                                            |                     |                  |                     |                  | AB                  | Fresenius Kabi        |                          |
| Amphoteric surfactants (amfotensiden)              | Tensodur 103                                | 2                    | 12                  |                                                       |                     |                                                                                                                            |                     |                  |                     |                  |                     | A                     | MFH Marienfelde          |
| Lye                                                | Lime-milk <sup>3</sup>                      |                      |                     |                                                       |                     |                                                                                                                            |                     | 20               | 6                   |                  |                     | A <sup>3</sup> B      |                          |

1 Not effective against myco-bacteria when service disinfecting, especially in the presence of blood.

2 Not suitable for disinfecting blood-contaminated or porous surfaces (e.g. raw wood).

3 Useless for tuberculosis; preparation of Lime-milk: 1 part dissolved lime (calcium hydroxide) + 3 parts water.

\* Checked for effectiveness on viruses in accordance with checking methods of the RKI (Federal Health Reporting 38 (1995) 242).

A: Suitable for killing vegetative bacterial germs including myco-bacteria as well as fungi, including fungal spores.

B: Suitable for deactivating viruses.

Table 4: Disinfectants

Kindly consult your specialist dealer if you have queries on matters related to disinfection; he will gladly assist you.

## 5 Guarantee

Excerpt from the "General Business Conditions":

(...)

5. The guarantee period for warranty claims is 24 months.

(...)

Claims are invalid where:

- the deviation from the agreed condition is unimportant
- usability is only minimally impaired
- structural changes to our products
- use of non-original parts
- wear and tear is normal
- faults are due to improper assembly or neglected maintenance
- faults are due to improper use
- errors are due to improper handling by customers or their contractual partners or service centres
- damage that arises after the transfer of risk as a result of faulty or negligent handling, excessive demands, unsuitable equipment, faulty assembly or that are due to particular external factors that are not assumed under the contract.

(...)

The terms of the guarantee may differ from country to country. Please contact your dealer for further information.

## 6 Disposal

The manufacturer is responsible for taking back and recycling the electric patient hoist and shall meet the requirements of European Directive 2002/96/EC on Waste Electrical and Electronic Equipment. You can find out from your local waste disposal company where you can take the electrical patient hoist to be recycled free of charge. It may not be disposed of with domestic waste.

Your dealer will also be able to advise you should you have any queries.

Packaging materials can be taken to disposal or recycling centers or to your specialist dealer.

## 7 Declaration of conformity

The patient hoist Eagle meet the requirements of European Directive  
- 93/42/EEC (Medical Products Directive)

and the product norms:  
- EN 10535: 2006

## 8 Maintenance plan

[illegible]

## 9 Disinfection book

[illegible]

**Abbreviations used in column 2 (reason):**

**V = Suspected infection    IF = Infection case**

**W = Repetition    I = Inspection**

## 10 Ordernumbers

### Eagle 620

|                                  |         |
|----------------------------------|---------|
| Casters (ø 100 mm)               | 1903403 |
| Braked casters (ø 100 mm)        | 1903402 |
| CBJ Home control box             | 1903418 |
| Adjustment motor LA 31C (6,000N) | 1903419 |
| Manual operation                 | 1903420 |
| Battery charger                  | 1903876 |
| Expander bracket                 | 16902xx |

### Eagle 625

|                                       |                    |
|---------------------------------------|--------------------|
| Casters (ø 100 mm)                    | 1903403            |
| Braked casters (ø 100 mm)             | 1903402            |
| E Control box                         | 1906137            |
| E Battery                             | 1906139            |
| E Construction plate                  | 1906230            |
| E Grey plug controlbox                | 1906232            |
| E Cable                               | 1906233            |
| E External charger CHJ2               | 1906423 (OPTIONAL) |
| E Construction plate external charger | 1906424 (OPTIONAL) |
| Adjustment motor LA 34 (10,000N)      | 1903887            |
| Manual operation                      | 1903420            |
| Expander bracket                      | 16902xx            |

### Slings

|                          |         |
|--------------------------|---------|
| Regular, Solid, Size M   | 1903824 |
| Regular, Solid, Size L   | 1903823 |
| Regular, Solid, Size XL  | 1905478 |
| Regular, Net, Size M     | 1905475 |
| Regular, Net, Size L     | 1905476 |
| Regular, Net, Size XL    | 1905477 |
| Comfort, Solid, Size M   | 1905359 |
| Comfort, Solid, Size L   | 1905251 |
| Comfort, Solid, Size XL  | 1905360 |
| Comfort, Net, Size M     | 1905479 |
| Comfort, Net, Size L     | 1905480 |
| Comfort, Net, Size XL    | 1905481 |
| Hygienic, Solid, Size M  | 1905357 |
| Hygienic, Solid, Size L  | 1905250 |
| Hygienic, Solid, Size XL | 1905358 |
| Hygienic, Net, Size M    | 1905482 |
| Hygienic, Net, Size L    | 1905483 |
| Hygienic, Net, Size XL   | 1905484 |

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ENGLISH

CONTRACTUAL WARRANTY

We offer 5 years of warranty on standard wheelchairs, lightweight wheelchairs 4 years. Electronic wheelchairs, tricycles, beds and other products: 2 years (batteries 6 months) and multiposition wheelchairs 3 years. This warranty is limited to the replacement of defective or spare parts.

APPLICATION CONDITIONS

In order to claim this warranty, part "B" of this card has to be given to your official Vermeiren dealer. The warranty is only valid when parts are replaced by Vermeiren in Belgium.

EXCEPTIONS

This warranty is not valid in case of:

- damage due to incorrect usage of the wheelchair,
- damage during transport,
- involvement in an accident,
- a dismount, modification or repair carried outside of our company and/or official Vermeiren dealership,
- normal wear of the wheelchair,
- non-return of the warranty card

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### **Vermeiren Nederland B.V.**